

WHILL Control System Protocol Specification

| | | |
|--|--|--------------|
|  | WHILL Control System Protocol Specification | Rev 9 |
|--|--|--------------|

Contents:

| | |
|---|----|
| 1. Preface | 2 |
| 2. Communication Overview..... | 2 |
| 3. Control Command | 3 |
| 3.1. Communication Format..... | 3 |
| 3.2. Control Commands | 4 |
| 3.2.1. StartSendingData command..... | 4 |
| 3.2.2. StopSendingData command..... | 5 |
| 3.2.3. SetPower | 5 |
| 3.2.4. SetJoystick command..... | 6 |
| 3.2.5. SetSpeedProfile command | 7 |
| 3.2.6. SetBatteryVoltageOut command (Model CR only)..... | 8 |
| 3.2.7. SetVelocity command | 8 |
| 4. WHILL state data..... | 10 |
| 4.1. Communication format..... | 10 |
| 4.2. Data set 0..... | 10 |
| 4.3. Data set 1..... | 11 |
| 5. Response Data..... | 13 |
| 6. Communication Electrical Specifications | 14 |
| 6.1. RS232C specification..... | 14 |
| 6.2. Minimum time between control commands..... | 15 |
| 6.3. Maximum interval of byte data in control command | 15 |
| 7. Sensor (Model CR only) | 16 |
| 7.1. 3D accelerometer/3D gyroscope sensor..... | 16 |
| 8. Connector..... | 16 |
| 8.1. Pin assignment | 16 |
| 8.1.1. Model CR..... | 16 |
| 8.1.2. Model CR2..... | 17 |
| 8.2. Battery voltage from connector (Model CR only) | 17 |
| 9. Appendix | 19 |

| | | |
|-------|---------------------|--------------|
| Rev 9 | Confidential | Page 1 of 21 |
|-------|---------------------|--------------|

| | |
|--|----|
| 9.1. Difference between Model CR and Model CR2 | 19 |
| 10. Revision History..... | 21 |

1. Preface

This document describes the specifications of the communication protocols and the control system of WHILL Model CR series provided by WHILL, Inc.

The Model CR series includes Model CR, Model CR2, Wheeled Robot Base, Electrical System Kit, and Omni-Platform. Wheeled Robot Base and Electrical System Kit have the same software specifications as Model CR2, and thus, all descriptions mentioned as Model CR2 in this document apply. Omni-Platform, a four-wheel-drive product, is out of scope of this document.




Figure 1 : WHILL Model CR series

2. Communication Overview

Communication interface between the host device and WHILL is RS232C.

In basic communication, the control command is sent from the host device to WHILL and WHILL takes actions according to the control command. The WHILL state data will be periodically sent from WHILL to the host device after WHILL receives the control command to send the WHILL state data. The WILL state data includes speed of WHILL, battery level,

battery current, 3D accelerometer/3D gyroscope sensor value and so on. Figure 1 shows the communication overview.




Figure 2 : Communication overview

3. Control Command

3.1. *Communication Format*

Figure 2 shows the communication format of the control command. The communication format includes protocol sign, data length, control command and checksum.

Protocol sign: At first, protocol sign is sent from the host device. WHILL recognizes the start of control command by protocol sign. Protocol sign is always 0xAF.

Data length: Data length shows byte size behind itself, not includes protocol sign and data length.

Control Command: Control command consists of command ID and command data.

Checksum: Checksum is the value of XOR of protocol sign, data length and control command (command ID and all command data).




Figure 3 : Communication format of the control command

3.2. Control Commands

The host device can control WHILL using control commands. Table 1 is the list of control command. Each command has a command ID and command data. The detail of each command is described from next section.

Table 1 : List of control commands

| Command | Commad ID | Description | Implementation | |
|-----------------------------|-----------|--|----------------|-----------|
| | | | Model CR | Model CR2 |
| StartSendingData | 0x00 | WHILL starts sending WHILL state data | ✓ | ✓ |
| StopSendingData | 0x01 | WHILL stops sending WHILL state data | ✓ | ✓ |
| SetPower | 0x02 | Power On/Off WHILL | ✓ | ✓ |
| SetJoystick | 0x03 | Disable/Enable joystick control and set joystick value | ✓ | ✓ |
| SetSpeedProfile | 0x04 | Set max speed, acceleration and deceleration | ✓ | ✓ |
| SetBatteryVoltageOut | 0x05 | Disable/Enable voltage from battery voltage pin of cable | ✓ | - |
| reserve | 0x06 | Don't use this Command ID | - | - |
| reserve | 0x07 | Don't use this Command ID | - | - |
| SetVelocity | 0x08 | Set velocity of WHILL | ✓ | ✓ |

3.2.1. StartSendingData command

StartSendigData command makes WHILL start sending the WHILL state data. Figure 3 shows StartSendingData command format.

Command ID: 0x00

D0: Data set number.

D0 = 0: Select data set 0

D1 = 1: Select data set 1

T0~T1: Interval of sending data [ms] (16bit value). T1:MSB8bit, T0:LSB8bit. Minimum interval is 10ms.

S0: Speed mode. Set speed mode in data set 0. See section 4.2 Data set 0 for more detail.

StartSendingData command




Figure 4 : StartSendingData format

3.2.2. StopSendingData command

StopSendigData command makes WHILL stop sending the WHILL state data. Figure 4 shows StopSendingData command format.

Command ID: 0x01

StopSendingData command




Figure 5 : StopSendingData format

3.2.3. SetPower

Function of SetPower command is to power ON/OFF WHILL. When WHILL finishes to power on after receiving SetPower command, WHILL sends Response data to the host device. Before receiving the response data, don't send any commands to WHILL. If the

response data isn't sent for 5ms, re-issue SetPower Command to WHILL. See section5 WHILL Response Data for more detail. Figure 5 shows SetPower command format.

Command ID: 0x02

P0: Select Power ON or OFF

P0 = 0: Power OFF

P0 = 1: Power ON

SetPower command




Figure 6 : SetPower command format

Note: After sending “POWER OFF”, the host device must wait for more than 5s to send “POWER ON”.

3.2.4. SetJoystick command

SetJoystick command has two functions. First function is to enable/disable the joystick control of a user who rides WHILL and disable/enable the joystick control of the host device. When the joystick control of a user is enabled (disabled), the joystick control of the host device is disabled (enabled). Second function is to set the value of joystick. The value is valid only when the joystick control of the host device is enabled. With a single SetJoystick command, the setting value is available for 200ms.Keep sending the command within 200ms if keeping WHILL moving more than 200ms.

Figure 6 shows SetJoystick command format.

Command ID: 0x03

U0: Disable/Enable joystick control of user/host device.

U0 = 0: Disable joystick control of user and enable joystick control of the host device

U0 = 1: Enable joystick control of user and disable joystick control of the host device

FB0 : Joystick value in front/back direction (-100~100)

LR0: Joystick value in left/right direction (-100~100)

SetJoystick command




Figure 7 : SetJoystick command format

3.2.5. SetSpeedProfile command

SetSpeedProfile command is to set max speed, acceleration and deceleration at forward movement, reverse movement and turn movement in each speed mode. Figure 7 shows SetSpeedProfile command format.

Command ID: 0x04

S1: Speed mode: Select speed mode to set

S1 = 0: Select speed mode when user uses joystick and “1” is displayed in indicator of WHILL

S1 = 1: Select speed mode when user uses joystick and “2” is displayed in indicator of WHILL

S1 = 2: Select speed mode when user uses joystick and “3” is displayed in indicator of WHILL

S1 = 3: Select speed mode when user uses joystick and “4” is displayed in indicator of WHILL

S1 = 4: Select speed mode when host device moves WHILL via RS232C

S1 = 5: Select speed mode when iPhone WHILL App moves WHILL

F_M1: Max speed at forward movement. The unit of value is 0.1km/h. If 30 is set, max speed is 3.0km/h. The range is (8~60).

F_A1: Acceleration at forward movement. The range is (10~90).

F_D1: Deceleration at forward movement. The range is (40~160).

R_M1: Max speed at reverse movement. The unit of value is 0.1km/h. If 30 is set, max speed is 3.0km/h. The range is (8~30).

R_A1: Acceleration at reverse movement. The range is (10~50).

R_D1: Deceleration at reverse movement. The range is (40~90).

T_M1: Max speed at turn movement. The unit of value is 0.1km/h. If 30 is set, max speed is 3.0km/h. The range is (8~35).

T_A1: Acceleration at turn movement. The range is (10~60).

T_D1: Deceleration at turn movement. The range is (40~160).

SetSpeedPrfile command




Figure 8 : SetForward command format

3.2.6. SetBatteryVoltageOut command (Model CR only)

SetBatteryVoltageOut command is to enable/disable battery voltage from battery voltage pin of cable connector. This voltage is connected to battery voltage. As default, battery voltage from battery voltage pin is enabled. See section 8.2 Battery voltage from connector about battery voltage pin.

Command ID: 0x05

V0: Enable/Disable voltage of battery voltage pin

V0 = 0: Disable

V0 = 1: Enable

SetBatteryVoltageOut command




Figure 9 : SetBatteryVoltageOut command format

3.2.7. SetVelocity command

SetVelocity command has two functions. First function is to enable/disable the joystick control of a user who rides WHILL and disable/enable the control of the host device. When the joystick control of a user is enabled (disabled), the control of the host device

is disabled (enabled). Second function is to set the value of WHILL velocity. The value is valid only when control of the host device is enabled. With a single SetVelocity command, the setting value is available for 200ms. Keep sending the command within 200ms if keeping WHILL moving more than 200ms.

[Difference from SetJoystick command]

Using SetVelocity command, the host device can send velocity value directly independently of speed profile setting. Also, a finer speed can be set using SetVelocity command than using SetJoystick command.

The acceleration is const when the host device uses SetVelocity command. The acceleration is fast (1.7m/s²) in order for WHILL to follow the target speed as quick as possible.

WHILL moves so quickly using SetVelocity command and so pay enough attention to use SetVelocity command. Basically, send this command to increase speed gradually.

Figure 9 shows SetVelocity command format.

Command ID: 0x06

U0: Disable/Enable joystick control of user/host device.

U0 = 0: Disable joystick control of user and enable joystick control of the host device

U0 = 1: Enable joystick control of user and disable joystick control of the host device

Y0~1: Velocity value in in front/back direction (16bit value). The range is (-500~1500) .The unit is 0.004[km/h]. Y1:MSB8bit, Y0:LSB8bit.

X0~1: Velocity value in in left/right direction (16bit value). The range is (-750~750) .The unit is 0.004[km/h]. X1:MSB8bit, X0:LSB8bit.

SetVelocity command




Figure 10 : SetVelocity command format

4. WHILL state data

4.1. Communication format

After receiving StartSendingData command, WHILL starts to send the WHILL state data to the host device. The host device selects the data set number of the WHILL state data. According to the data set number, contents which WHILL sends changes. The WHILL state data conforms to the communication format shown in Figure 10.

Protocol sign: At first, protocol sign is sent from WHILL. Protocol sign indicates the start of the WHILL state data. Protocol sign is always 0xAF.

Data length: Data length shows byte size behind itself, not includes protocol sign and data length.

WHILL state data: WHILL state data consists of data set number and information.

Checksum: Checksum is the value of XOR of protocol sign, data length and WHILL state data.




Figure 11 : Communication format of WHILL state data

4.2. Data set 0

- 1.1.1 When the host device select data set 0 in StartSendingData command, values of Table 2 will be sent from WHILL. This data is about speed profile. This data set will be changed only when the host device modifies the setting of WHILL. So, it is better to read the data set 0 once when WHILL power-on or the host device changes the settings of max speed, acceleration and deceleration. See section 3.2.5 SetSpeedProfile command to know items in table 2.

Table 2 : Data set 0 contents

| Information Number | Value (8bit) | Description | Implementation | |
|--------------------|--------------|-------------|----------------|-----------|
| | | | Model CR | Model CR2 |

| | | | | |
|---|-------------------|----------------------------------|---|---|
| 0 | SPEED_MODE | Speed mode to send | ✓ | ✓ |
| 1 | FORWARD_SPEED_MAX | Max speed at forward movement | ✓ | ✓ |
| 2 | FORWARD_ACCEL | Acceleration at forward movement | ✓ | ✓ |
| 3 | FORWARD_DECEL | Deceleration at forward movement | ✓ | ✓ |
| 4 | REVERSE_SPEED_MAX | Max speed at reverse movement | ✓ | ✓ |
| 5 | REVERSE_ACCEL | Acceleration at reverse movement | ✓ | ✓ |
| 6 | REVERSE_DECEL | Deceleration at reverse movement | ✓ | ✓ |
| 7 | TURN_SPEED_MAX | Max speed at turn movement | ✓ | ✓ |
| 8 | TURN_ACCEL | Acceleration at turn movement | ✓ | ✓ |
| 9 | TURN_DECEL | Deceleration at turn movement | ✓ | ✓ |

4.3. Data set 1

When the host device selects data set 1 in StartSendingData command, values of Table 3 will be sent from WHILL. This data set is changing continuously such as acceleration sensor value, battery current and so on.

Table 3 : Data set 1 contents

| Information Number | Value (8bit) | Description | Implementation | |
|--------------------|------------------|---|----------------|-----------|
| | | | Model CR | Model CR2 |
| 0 | ACC_X (MSB 8bit) | Acceleration X axis (MSB 8bit) See Note1 | ✓ | - |
| 1 | ACC_X (LSB 8bit) | Acceleration X axis (LSB 8bit) See Note1 | ✓ | - |
| 2 | ACC_Y (MSB 8bit) | Acceleration Y axis (MSB 8bit) See Note1 | ✓ | - |
| 3 | ACC_Y (LSB 8bit) | Acceleration Y axis (LSB 8bit) See Note1 | ✓ | - |
| 4 | ACC_Z (MSB 8bit) | Acceleration Z axis (MSB 8bit) See Note1 | ✓ | - |
| 5 | ACC_Z (LSB 8bit) | Acceleration Z axis (LSB 8bit) See Note1 | ✓ | - |
| 6 | GYR_X (MSB 8bit) | Angular rate X axis (MSB 8bit) See Note2 | ✓ | - |
| 7 | GYR_X (LSB 8bit) | Angular rate X axis (LSB 8bit) See Note2 | ✓ | - |
| 8 | GYR_Y (MSB 8bit) | Angular rate Y axis (MSB 8bit) See Note2 | ✓ | - |
| 9 | GYR_Y (LSB 8bit) | Angular rate Y axis (LSB 8bit) See Note2 | ✓ | - |

| | | | | |
|----|---------------------------------|---|---|---|
| 10 | GYR_Z (MSB 8bit) | Angular rate Z axis (MSB 8bit) See Note2 | ✓ | - |
| 11 | GYR_Z (LSB 8bit) | Angular rate Z axis (LSB 8bit) See Note2 | ✓ | - |
| 12 | JOY_FRONT | Joystick value in front/back direction (-100 ~ +100). Value of real joystick, not value set by the host device | ✓ | - |
| 13 | JOY_SIDE | Joystick value in left/right direction (-100 ~ +100). Value of real joystick, not value set by the host device | ✓ | - |
| 14 | BATTERY_POWER | Battery level. 0 ~ 100% | ✓ | ✓ |
| 15 | BATTERY_CURRENT (MSB 8bit) | Battery current (MSB 8bit) See Note3 | ✓ | ✓ |
| 16 | BATTERY_CURRENT (LSB 8bit) | Battery current (LSB 8bit) See Note3 | ✓ | ✓ |
| 17 | RIGHT_MOTOR_ANGLE (MSB 8bit) | Angle of right motor (MSB 8bit). See Note4 | ✓ | ✓ |
| 18 | RIGHT_MOTOR_ANGLE (LSB 8bit) | Angle of right motor (LSB 8bit). See Note4 | ✓ | ✓ |
| 19 | LEFT_MOTOR_ANGLE (MSB 8bit) | Angle of left motor (MSB 8bit). See Note4 | ✓ | ✓ |
| 20 | LEFT_MOTOR_ANGLE (LSB 8bit) | Angle of left motor (LSB 8bit). See Note4 | ✓ | ✓ |
| 21 | RIGHT_MOTOR_SPEED (MSB 8bit) | Right motor speed value (MSB 8bit). See Note5 | ✓ | ✓ |
| 22 | RIGHT_MOTOR_SPEED (LSB 8bit) | Right motor speed value (LSB 8bit). See Note5 | ✓ | ✓ |
| 23 | LEFT_MOTOR_SPEED (MSB 8bit) | Left motor speed value (MSB 8bit). See Note5 | ✓ | ✓ |
| 24 | LEFT_MOTOR_SPEED (LSB 8bit) | Left motor speed value (LSB 8bit). See Note5 | ✓ | ✓ |
| 25 | POWER_ON | State of WHILL power. ON:1, OFF:0 | ✓ | ✓ |
| 26 | SPEED_MODE_INDICATOR | Speed mode displayed on indicator of WHILL | ✓ | ✓ |
| 27 | ERROR | Error code on WHILL. If there is no error, '0' is sent. | ✓ | ✓ |
| 28 | ANGLE_DETECT_COUNTER | Counter of timing to detect motor angle. | ✓ | ✓ |

Note1: Acceleration value must be calculated below in order to be converted into unit [mg].

$$V_{acc} \times 0.122 = a$$

V_{acc} : Value about acceleration which is reported from WHILL. 16bit signed value.

a : Acceleration value which unit is [mg].

Note2: Angular rate value must be calculated below in order to be converted into unit [mdps].

| | | |
|--|--|--------------|
|  | WHILL Control System Protocol Specification | Rev 9 |
|--|--|--------------|

$$V_{ang} \times 4.375 = \omega$$

V_{ang} : Value about angular rate which is reported from WHILL. 16bit signed value.
 ω : Angular rate value which unit is [mdps].

Note3: BATTERY_CURRENT is 16bit signed value. The unit is 2mA.
Example: 0x035 represents 106mA; 0xFF97 represents -210mA.
The value will be 0 if battery current is in ± 75 mA.
Sampling rate of battery current is 4Hz.

Note4: RIGHT_MOTOR_ANGLE and LEFT_MOTOR_ANGLE is 16bit signed value. The unit is 0.001rad. The range is $\pm \pi$ [rad].
Example: 0x0600 represents 1.536rad

Note5: The value of the item that Implementation is “-” is always 0x00.

Note5: RIGHT_MOTOR_SPEED and LEFT_MOTOR_SPEED is 16bit signed value. The unit is 0.004km/h.
Example: 0x01F4 represents 2km/h.
If higher accuracy of motor speed is needed, calculate motor speed from RIGHT_MOTOR_ANGLE(LEFT_MOTOR_ANGLE) and ANGLE_DETECT_COUNTER instead of using the value of RIGHT_MOTOR_SPEED(LEFT_MOTOR_SPEED).

Note6: ANGLE_DETECT_COUNTER indicates a timing of detecting motor angles. The unit is [ms].
The max of this value is 200. After this value reaches 200, this value will be 0 in next timing of detecting motor angles.
This value is used in order to improve the accuracy of odometry.
For Example: When the following data is received sequentially, right motor rotates 0.088rad in 11ms.

RIGHT_MOTOR_ANGLE: 1.536rad and ANGLE_DETECT_COUNTER: 110
RIGHT_MOTOR_ANGLE: 1.624rad and ANGLE_DETECT_COUNTER: 121

5. Response Data

When WHILL finishes to powers on after receiving SetPower command, WHILL sends the response data to the host device. The response data is always 0x52. Figure 11 shows the response data format.

Protocol Sign: At first, protocol sign is sent from WHILL. Protocol sign is always 0xAF.

Data length: it shows byte size behind itself, not includes protocol sign and data length. It's '2' at WHILL response data.

| | | |
|-------|--------------|---------------|
| Rev 9 | Confidential | Page 13 of 21 |
|-------|--------------|---------------|

Response: It's always 0x52

Checksum: it is the value of XOR of protocol sign, data length and response.

WHILL response data




Figure 12 : WHILL response data format

After sending SetPower command to power on WHILL, don't send any commands to WHILL before receiving the response data. If the response data isn't sent for 15ms, re-issue SetPower command to WHILL (Figure 12).




Figure 13 : Timing of Response data

6. Communication Electrical Specifications

6.1. RS232C specification

The interface of WHILL's CPU is RS232C base on Table 4. Figure 13 shows RS232C diagram.

Table 4 : UART specification

| | Model CR | Model CR2 |
|-------------|----------|-----------|
| Baud rate | 38400 | 38400 |
| Parity bit | No | No |
| Data length | 8bit | 8bit |
| Stop bit | 2bit | 2bit |



Figure 14 : RS232C diagram

6.2. *Minimum time between control commands*

The host device must wait for 2ms to the send next control command after sending the previous control command (Figure 14).




Figure 15 : Minimum time between control commands

6.3. *Maximum interval of byte data in control command*

The interval of byte data which constitutes a single control command must be less than 5ms. If the interval is more than 5ms, the previous 1 byte data will be ignored.




Figure 16 : Maximum interval of byte data in control command

| | | |
|--|--|--------------|
|  | WHILL Control System Protocol Specification | Rev 9 |
|--|--|--------------|

7. Sensor (Model CR only)

7.1. 3D accelerometer/3D gyroscope sensor

3D accelerometer/3D gyroscope sensor is included in WHILL. Acceleration sampling rate is 104Hz and acceleration range is $\pm 4g$. Angular rate sampling rate is 104Hz and angular rate range is $\pm 125\text{dps}$. Figure 16 shows direction of the detectable accelerations and angular rate.

For the sensor, LSM6DS3 is included in WHILL. See LSM6DS3 datasheet for more detail.



Figure 17 : Direction of the detectable accelerations and angular rate

8. Connector


8.1. Pin assignment

8.1.1. Model CR

08R-JWPF-VSLE-D (JST) is used for interface connector. Figure 17 shows pin assignment.

| | | |
|-------|--------------|---------------|
| Rev 9 | Confidential | Page 16 of 21 |
|-------|--------------|---------------|

| | |
|-----|------|
| TXD | RXD |
| GND | BATT |
| GND | BATT |
| GND | BATT |



TXD : TX pin. Data from WHILL
 RXD: Rx pin. Data to WHILL
 GND: GND
 BATT: Battery voltage pin. This pin connected with battery voltage.

Figure 18 : Pin assignment of connector

8.1.2. Model CR2

JEZ-9S-3 (JST) is used for interface connector. Figure 18 shows pin assignment.




Figure 19 : Pin assignment of connector

8.2. **Battery voltage from connector (Model CR only)**

BATT pin in Figure 17 is connected to battery voltage. Battery voltage range from about 19V ~ 30V according to remaining battery level, battery current and battery temperature. This voltage output can be enable / disable by SetBatteryVoltageOut command.

The current from this BATT pin must be under 3A. The fuse between BATT pin and battery will blows if over current flows.

9. Appendix

9.1. *Difference between Model CR and Model CR2*

Below, check-mark ("✓") means IMPLEMENTED.

Table 5 : List of External Pins

| External Pins | Model CR | Model CR2 |
|---------------------------------|----------|-----------|
| RS232C (TX, RX) | ✓ | ✓ |
| Battery voltage output (24Vout) | ✓ | - |

Table 6 : List of UART specifications

| UART specifications | Model CR | Model CR2 |
|---------------------|----------|-----------|
| Baud rate | 38400 | 38400 |
| Parity bit | No | No |
| Data length | 8bit | 8bit |
| Stop bit | 2bit | 2bit |

Table 7 : List of Control Commands

| Command ID | Command | Model CR | Model CR2 |
|------------|----------------------|----------|-----------|
| 0 | StartSendingData | ✓ | ✓ |
| 1 | StopSendingData | ✓ | ✓ |
| 2 | SetPower | ✓ | ✓ |
| 3 | SetJoystick | ✓ | ✓ |
| 4 | SetSpeedProfile | ✓ | ✓ |
| 5 | SetBatteryVoltageOut | ✓ | - |
| 8 | SetVelocity | ✓ | ✓ |

Table 8 : List of WHILL State Data (Data Set 0)

| Information Number | Value (8bit) | Model CR | Model CR2 |
|--------------------|-------------------|----------|-----------|
| 0 | SPEED_MODE | ✓ | ✓ |
| 1 | FRONT_SPEED_MAX | ✓ | ✓ |
| 2 | FRONT_ACCEL | ✓ | ✓ |
| 3 | FRONT_DECEL | ✓ | ✓ |
| 4 | REVERSE_SPEED_MAX | ✓ | ✓ |
| 5 | REVERSE_ACCEL | ✓ | ✓ |
| 6 | REVERSE_DECEL | ✓ | ✓ |
| 7 | TURN_SPEED_MAX | ✓ | ✓ |

| | | | |
|---|------------|---|---|
| 8 | TURN_ACCEL | ✓ | ✓ |
| 9 | TURN_DECEL | ✓ | ✓ |

Table 9 : List of WHILL State Data (Data Set 1)

| Information Number | Value (8bit) | Model CR | Model CR2 |
|--------------------|------------------------------|----------|-----------|
| 0 | ACC_X (MSB 8bit) | ✓ | - |
| 1 | ACC_X (LSB 8bit) | ✓ | - |
| 2 | ACC_Y (MSB 8bit) | ✓ | - |
| 3 | ACC_Y (LSB 8bit) | ✓ | - |
| 4 | ACC_Z (MSB 8bit) | ✓ | - |
| 5 | ACC_Z (LSB 8bit) | ✓ | - |
| 6 | GYR_X (MSB 8bit) | ✓ | - |
| 7 | GYR_X (LSB 8bit) | ✓ | - |
| 8 | GYR_Y (MSB 8bit) | ✓ | - |
| 9 | GYR_Y (LSB 8bit) | ✓ | - |
| 10 | GYR_Z (MSB 8bit) | ✓ | - |
| 11 | GYR_Z (LSB 8bit) | ✓ | - |
| 12 | JOY_FRONT | ✓ | - |
| 13 | JOY_SIDE | ✓ | - |
| 14 | BATTERY_POWER | ✓ | ✓ |
| 15 | BATTERY_CURRENT (MSB 8bit) | ✓ | ✓ |
| 16 | BATTERY_CURRENT (LSB 8bit) | ✓ | ✓ |
| 17 | RIGHT_MOTOR_ANGLE (MSB 8bit) | ✓ | ✓ |
| 18 | RIGHT_MOTOR_ANGLE (LSB 8bit) | ✓ | ✓ |
| 19 | LEFT_MOTOR_ANGLE (MSB 8bit) | ✓ | ✓ |
| 20 | LEFT_MOTOR_ANGLE (LSB 8bit) | ✓ | ✓ |
| 21 | RIGHT_MOTOR_SPEED | ✓ | ✓ |
| 22 | RIGHT_MOTOR_SPEED | ✓ | ✓ |
| 23 | LEFT_MOTOR_SPEED | ✓ | ✓ |
| 24 | LEFT_MOTOR_SPEED | ✓ | ✓ |
| 25 | POWER_ON | ✓ | ✓ |
| 26 | SPEED_MODE_INDICATOR | ✓ | ✓ |
| 27 | ERROR | ✓ | ✓ |
| 28 | ANGLE_DETECT_COUNTER | ✓ | ✓ |

10. Revision History

| Revision | Date | Note |
|----------|------------|--|
| Rev 9 | 2024/03/27 | Add description of “WHILL Model CR series” in Preface. |
| Rev 8 | 2023/11/10 | Add description of Model CR2 (Table1, 2, 3, 4). Add Pin assignment of Model CR2. Add Appendix (Difference between Model CR and Model CR2). Modify waiting time of response. |
| Rev 7 | 2019/12/25 | Add ANGLE_DETECT_COUNTER in Data set 1. Modify Table3:Data set 1 contents. Add SETVelocity command. |
| Rev 6 | 2019/1/10 | Modify sentence. Modify description about SetPower and Response data according to new FW. |
| Rev 5 | 2018/09/14 | Modify Figure 15 |
| Rev 4 | 2018/01/31 | Delete EmergencyBrake command. Modify range of max speed, acceleration and deceleration in setSpeedProfile command. |
| Rev 3 | 2017/10/30 | Add SPEED_MODE_INDICATOR in Data set 1. Add EmergencyBrake command. |
| Rev 2 | 2017/08/23 | Modify figures of commands. Modify sentence. |
| Rev 1 | 2017/08/22 | |